



Step 1: milling

BioKey Core Size: 36 x 47 x 29 mm (width x height x depth). Milling Size: >37x48x35 mm, **just one Milling** is necessary, construction height: > 120cm, stainless steel panel: 44,6x75,4 x 1.5mm. Screw (2x, Type: C2.9x22H, distance: 61 mm).





Step 2: cabling

Cabeling in the door, to use just for manufacturer 8 ~ 24 V DC, stabilized potentialfree Contact, max. 15mA Power Supply 8-24V DC + Yellow BioKey-Inside Modul brown + **BioKey Inside** White Motor Lock 1) Green Potentialfree Contact Motor Lock White to + (Plus) KFV GU A-O + (Plus) **Cabeling with Motor Lock** Fuhr - (Minus)

Step 3: Autotest

To test the functionalities (without Fingerprints) with the IR-Keypad: **0 OK**. If OK, green/red LEDs glow constantly (after about 10s).



IDENCOM Germany GmbH, Schillerstr. 56-58, 10627 Berlin, Germany Tel. +49-30-398839 310, Fax: +49-30-398839 318 E-mail: <u>support@idencom.com</u> IDENCOM Germany GmbH- March 15 – Printed in Germany- All rights reserved BioKey® is a registered trademark of IDENCOM AG (Schweiz / Switzerland). All other names are the property of there respective owners.